

# FreeStyle LibreLink App Setup Overview



Download FreeStyle LibreLink from App Store



**YOU CAN DO IT**  
WITHOUT LANCETS<sup>1</sup>

Android



iOS



## 1 Assemble and Apply Sensor to the back of your upper arm

### Get Ready

#### STEP 1



Select site on back of upper arm. Do not use other sites as these are not approved and may result in inaccurate glucose readings.

**Note:** Avoid scars, moles, stretch marks, lumps, and insulin injection sites. To prevent skin irritation, rotate sites between applications.

#### STEP 2



Clean site with alcohol wipe. Allow site to dry before proceeding.

#### STEP 3



Peel lid completely off Sensor Pack. Unscrew cap from Sensor Applicator.

**CAUTION:** Sensor codes must match on Sensor Pack and Sensor Applicator or glucose readings will be incorrect.



### Prepare Sensor Applicator

#### STEP 4



Line up dark mark on Sensor Applicator with dark mark on Sensor Pack. On a hard surface, press down firmly on Sensor Applicator until it comes to a stop.

#### STEP 5



Lift Sensor Applicator out of Sensor Pack.

#### STEP 6



Sensor Applicator is ready to apply Sensor.

**CAUTION:** Sensor Applicator now contains a needle. Do not touch inside Sensor Applicator or put it back into Sensor Pack.

### Apply Sensor


#### STEP 7



Place Sensor Applicator over site and push down firmly to apply Sensor.

**CAUTION:** Do not push down on Sensor Applicator until placed over prepared site to prevent unintended results or injury.

#### STEP 8



Gently pull Sensor Applicator away from your body.

#### STEP 9



Make sure Sensor is secure. Discard used Sensor Applicator and Sensor Pack according to local regulations.


## 2 Set up App

### STEP 1



Check your smartphone is connected to a WiFi or cellular network. Tap **Sign In** or **GET STARTED NOW**. Review the legal information and follow instructions to set up a LibreView account or login to your existing account.

### STEP 2



Enter your personal settings. Tap NEXT after each screen.

### STEP 3



Review important information about glucose readings, when to do a blood glucose test, and where to apply a Sensor. Tap NEXT to proceed until you see the Scan New Sensor screen.

## 3 Use the App to scan your Sensor

### How to Scan

Tap scan button (📶). You can either tap the blue box on the Home Screen or (📶) at the top right. Hold top of iPhone near Sensor until you get a tone and/or vibration.


### STEP 1



Scan your Sensor to start it.


**Note:** If you want to use both the Reader and the App to check your glucose, you must start the Sensor with the Reader first.

### STEP 2




The Sensor can be used to check your glucose after the start-up period. While the Sensor is starting up, you can navigate away from the App.

### STEP 3



After the start-up period, scan your Sensor to get your glucose reading.

## Reports of your glucose history



- Tap ☰ at the top left of the screen to access your reports.
- Tap 📷 to share a screenshot of the report.
- Tap ⓘ to view a description of the report.



1. Scanning the sensor does not require lancets.  
 The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink requires registration with LibreView.  
 Simulated data for illustrative purposes only. Not real patient or data.  
 Images are for illustrative purposes only. Not real patient or Health Care Professional.  
 © 2021 Abbott  
 FreeStyle, Libre, and related brand marks are marks of Abbott.  
 ADC-18402 V2.0 5/21 ZHF3917 - 10/05/2022