

# Get connected and start sharing today!

Now it's even easier to share your glucose data with your doctor



**YOU CAN DO IT REMOTELY**




## CONNECT TO A PRACTICE using the FreeStyle LibreLink app<sup>1</sup>

Connect to a Practice is a FreeStyle LibreLink feature that lets you easily share your glucose data with your doctor's practice through LibreView<sup>2</sup>



FreeStyle LibreLink



## Discover glucose patterns and trends

Your glucose data is automatically uploaded to LibreView, a cloud-based system that enables your doctor to access reports, on demand



LibreView



Images and simulated data for illustrative purposes only. Not real patient or data.



Scan the QR code with your iPhone camera or tap your Android device to download the FreeStyle LibreLink app<sup>1</sup>



## How to connect to your Doctor's clinic remotely

In the FreeStyle LibreLink app<sup>1</sup>, follow the steps below, and enter your doctor's LibreView<sup>2</sup> Practice ID.

- 1** Click **Connected Apps**
- 2** Select **LibreView**
- 3** Click **Connect to a Practice**
- 4** Enter **Practice ID** shared by your Doctor

Practice ID:



Disclaimers:

1- The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView, a service provided by Abbott and Newyu, Inc.  
2- The LibreView website is only compatible with certain operating systems and browsers. Please check [www.LibreView.com](http://www.LibreView.com) for additional information. The user must be connected to data services for glucose data to automatically upload to LibreView.

©2021 Abbott. FreeStyle, Libre, and related brand marks are marks of Abbott.

ADC-32192 v2.0 01/2021  
MOHAP - ATGA1FEY - Exp 25/4/22