



# SET UP YOUR PRACTICE TO MONITOR CONNECTED PATIENTS

87% of patients sharing results with their physicians stay active on FreeStyle Libre 14 day sensors after 12 months.\*

## NOW YOU KNOW



  
**FreeStyle Libre 14 day**  
FLASH GLUCOSE MONITORING SYSTEM

  
*life. to the fullest.®*  
**Abbott**

The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.

\* In comparison to FreeStyle LibreLink app patients not sharing with their physicians. Data on file.

# Start with LibreView<sup>†</sup>

## LibreView

LibreView is a secure,<sup>‡</sup> cloud-based system that enables more effective diabetes management by giving practices visibility to their patients' FreeStyle Libre 14 day system glucose data.



### Easy upload

Glucose data are automatically<sup>‡</sup> uploaded from a patient's FreeStyle LibreLink app<sup>§</sup> to LibreView and can be manually uploaded from the FreeStyle Libre 14 day reader.



### Remote access

Access patient data securely and on demand by logging into any internet-connected computer. (See LibreView.com for system requirements.)



### Actionable insights

View trends, patterns, and easy-to-read standardized reports, including Time in Range, to make more informed treatment decisions.



Get started at [www.LibreView.com](http://www.LibreView.com)

## Components of monitoring your connected patients

### Enabling your practice to connect

It starts with creating a LibreView account and Practice ID.

### Inviting your FreeStyle Libre 14 day system patients to share their glucose data with you

Data sharing populates insightful glucose reports.

### Discovering actionable insights from shared data

Using LibreView glucose reports can help you to:

- Observe your patients' progress remotely
- Have more productive conversations with your patients
- Make more informed treatment decisions

<sup>†</sup> The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

<sup>‡</sup> The user's device must have internet connectivity for glucose data to automatically upload to LibreView and to transfer to connected LibreLinkUp app users.

<sup>§</sup> The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.


<sup>¶</sup> LibreView is ISO27001/27018/27701 certified and HITRUST CSF Certified.

# Set Up Your LibreView Practice

## + Create your LibreView Practice ID

Creating a LibreView Practice generates a LibreView Practice number that patients can use to share glucose data from anywhere.† Follow these steps to get started:

- 1 Create or sign into your LibreView Practice account at [www.LibreView.com](http://www.LibreView.com)
- 2 Select the  LibreView menu icon
- 3 Select **Create a New Practice**, then select **Continue**
- 4 Complete the required information, then select **Create a New Practice** to create your LibreView Practice and generate your LibreView Practice ID.




## ✎ Customize your LibreView Practice ID

Now that you've created a LibreView Practice ID, you can customize it, if desired, to make it easier to remember and share with your patients.

To customize your LibreView Practice ID:

- 1 Select the  LibreView menu icon
- 2 Select the **My Practices** drop-down list
- 3 Select the LibreView Practice you would like to view
- 4 Select **Edit** to change your LibreView Practice settings or LibreView Practice ID




† The user's device must have internet connectivity for glucose data to automatically upload to LibreView and to transfer to connected LibreLinkUp app users.  
ADC-20721 v4.0 2/22

# OPTION 1: Invite patients to share glucose data with Practice ID

## Provide your LibreView Practice ID to patients

Once you've created a Practice in LibreView,<sup>†</sup> you can invite your patients using the FreeStyle Libre 14 day system to share their glucose data remotely. At any time, refer back to your LibreView Practice ID by selecting  **LibreView**. Your LibreView Practice ID will appear under **My Practices**, below your practice name.



LibreView

Patients can enable data sharing with your LibreView Practice ID using the FreeStyle LibreLink app.<sup>§, #</sup>

Your Practice ID is available in the  **LibreView** menu under **Practices**.

This method can help patients initiate a connection themselves if receiving an email or app invitation isn't preferred.



Patients can enter your LibreView Practice ID by:

- 1 Selecting the app's  menu icon
- 2 Tapping **Connected Apps**
- 3 Tapping **LibreView**
- 4 Tapping **Connect to a Practice**
- 5 Entering your LibreView Practice ID and following the prompts to connect



Patients using FreeStyle Libre 14 day readers can enter your LibreView Practice ID by:

- 1 Signing into or creating a patient LibreView account
- 2 Navigating to **Account Settings** from the  **LibreView** menu icon
- 3 Selecting **My Practices**
- 4 Entering your LibreView Practice ID and clicking **Add**

They can then plug in their reader and follow the prompts from LibreView to download the LibreView device drivers and upload glucose data.

**Need help?**  
Visit the **Support** section of [www.LibreView.com](http://www.LibreView.com)

<sup>†</sup> The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

<sup>§</sup> The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.

<sup>#</sup> The FreeStyle LibreLink app and the FreeStyle Libre 14 Day reader have similar but not identical features. Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

# OPTION 2: Invite patients to share glucose data with their email




## Invite patients using their email



### LibreView

Invite patients to connect from LibreView†

- 1 Search for an existing patient email address, or click the icon at the top of the navigator bar to add a new one to your LibreView Practice
- 2 Select the patient's profile, then select **Invite** and use the email address associated with their account
- 3 Your invitation will appear inside the patient's FreeStyle LibreLink app§,¶ and in their email inbox



Patients can accept your invitation directly in the FreeStyle LibreLink app by:

- 1 Selecting the app's menu icon
- 2 Tapping **Connected Apps**
- 3 Tapping **LibreView**
- 4 Completing the prompts to connect



Patients using FreeStyle Libre 14 day readers can accept your invitation by:

- 1 Receiving your invitation to connect in their email inbox
- 2 Clicking **Accept Invitation** in the invitation
- 3 Following prompts to sign in or create a patient LibreView account

Afterward, they can plug in their FreeStyle Libre 14 day reader and follow the prompts from LibreView to download the LibreView device drivers and upload glucose data.

**Need help?**  
Visit the **Support** section of [www.LibreView.com](http://www.LibreView.com)

† The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

§ The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.

¶ The FreeStyle LibreLink app and the FreeStyle Libre 14 Day reader have similar but not identical features. Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

# Discover actionable insights in detailed glucose reports




## Gather insights

After you establish a data sharing connection with your FreeStyle Libre 14 day system patients, you can access insightful glucose reports for each, enabling more productive conversations and more informed treatment decisions. If patients are using both devices,† data will be seamlessly merged into reports.



## LibreView

Standardized reports:



† The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

‡ The user's device must have internet connectivity for glucose data to automatically upload to LibreView and to transfer to connected LibreLinkUp app users.

§ The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.

# The FreeStyle LibreLink app and the FreeStyle Libre 14 Day reader have similar but not identical features. Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

# Recommend the FreeStyle LibreLink app<sup>s,#</sup> and easily connect with your patients through LibreView<sup>†</sup>



LibreView

Healthcare professional support

Call 1-754-210-8584, ext 3, or visit the Support section on LibreView.com.



Patient support

Contact our Customer Care Team at 1-855-632-8658  
7 days a week, 8 AM to 8 PM ET, excluding holidays



† The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

§ The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check our website for more information about device compatibility before using the app. Use of the FreeStyle LibreLink app requires registration with LibreView.

# The FreeStyle LibreLink app and the FreeStyle Libre 14 Day reader have similar but not identical features. Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

Reference: 1. Data on file. Abbott Diabetes Care.

#### Important Safety Information

Failure to use FreeStyle Libre 14 day system as instructed in labeling may result in missing a severe low or high glucose event and/or making a treatment decision, resulting in injury. If readings do not match symptoms or expectations, use a fingerstick value from a blood glucose meter for treatment decisions. Seek medical attention when appropriate or contact Abbott at 855-632-8658 or [FreeStyleLibre.us](http://FreeStyleLibre.us) for safety info.

The circular shape of the sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners.

©2022 Abbott. All rights reserved. ADC-20721 v4.0 2/22