



FreeStyle

Libre



LibreView

Easy insights¹.

Use LibreView² to share glucose readings with your healthcare team for more effective consultations.

Learn more at **FreeStyle.Abbott**

Images are for illustrative purposes only. Not actual patient.

1. Unger, J. *Postgrad Med.* (2020): <https://doi.org/10.1080/00325481.2020.1744393>. 2. The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose device data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.




LibreView

Clear, intuitive reports for you and your healthcare professional

SETTING UP YOUR ACCOUNT

- Step 1 Click on **sign up** and select **LibreView Account** for patients
- Step 2 Enter the requested **Account information**
- Step 3 Verify your email
- Step 4 Complete 2-factor authentication



Sign up at www.LibreView.com

UPLOADING YOUR DATA



Once you have installed the LibreView device drivers, you can upload data into your LibreView account **at home, anytime.**

- Step 1 Connect your Libre 2 or Libre 3 reader to a PC via your yellow cable
- Step 2 Press to **Begin Upload**
- Step 3 Follow the on-screen instructions to finish installation
- Step 4 To ensure all data is captured, upload your reader at least once every 90 days

Note: If you are using a compatible mobile app within the FreeStyle Libre portfolio to monitor your glucose, your data is automatically uploaded to LibreView when an internet connection is available.


Images are for illustrative purposes only. Not actual patient data.


The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose device data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice.

LibreView

VIEWING REPORTS¹

Once you have uploaded your data, you can view your reports as follows:

- Click on 
- Select **Glucose reports**
- Choose the report(s) you wish to view from the left-hand menu




CONNECTING TO YOUR PRACTICE

You can link your LibreView account with your healthcare team to easily share your data and optimise appointments.

OPTION 1

- Ask your healthcare professional to invite you to share data with them
- You will receive an email invitation to join their practice
- Click on **Accept**

OPTION 2

- Ask your healthcare professional for their unique LibreView Practice ID
- Log into your LibreView account
- Select  in the top right of the screen
- Click on **Account Settings**, then **My Practices**
- Enter the Practice ID and click **Add**

Images are for illustrative purposes only. Not actual patient or data.

¹ The LibreView website is only compatible with certain operating systems and browsers. Please check www.LibreView.com for additional information.



LibreView – a secure¹, cloud-based diabetes management system

See and share your complete glycaemic profile²



Clear, intuitive reports

Identify glucose patterns and trends thanks to eight easy-to-read reports so you can make informed decisions about your health.



Easy, remote access

View your glucose data online anytime, anywhere, from any internet-connected device.



Connected care

Easily share your reports with your healthcare team, even between appointments, to enable more productive consultations.

1. LibreView is ISO27001/27018/27701 certified and HITRUST CSF certified. 2. For a complete glycaemic picture, scan once every eight hours during signal loss when using the FreeStyle LibreLink app. Reader users need to scan their sensors once every eight hours for a complete glycaemic picture.

© 2025 Abbott. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott.
ADC-87728 v2.0 03/25.

