

The FreeStyle Libre 14 day system is an accurate,<sup>1</sup> convenient, and easy<sup>2</sup> way to get glucose readings and trends.

## NOW YOU KNOW



### Two components:

#### Innovative sensor technology



- Start taking glucose readings just one hour after sensor application<sup>\*</sup>
- Small, discreet,<sup>2</sup> and water resistant<sup>†</sup>
- Can be worn for up to 14 days
- No fingerstick calibration required<sup>§</sup>
- Sensor automatically measures glucose readings day and night<sup>||</sup>



#### Compact, lightweight reader

- Painless,<sup>2</sup> one-second scan—even through clothing<sup>¶</sup>
- Stores up to 90 days of glucose data
- Backlit touchscreen for visibility

### How to get started with the FreeStyle Libre 14 day system

#### 1 Get a prescription

Request a written prescription of 1 reader and 2 sensors (a 28-day supply) with 12 refills (a 1-year supply).

#### 2 Go to a pharmacy

Fill your prescription at a participating pharmacy.<sup>\*\*</sup> Your pharmacy may have to order the product for you.

The FreeStyle Libre 14 day system is available with a prescription at participating pharmacies.<sup>\*\*</sup>

\* Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

† Glucose readings are not available during 1-hour warmup, when sensor is too hot or too cold or when you see an error or "LO" or "HI" message or no current glucose reading.

‡ Sensor is water resistant in up to 1 meter (3 feet) of water. Do not immerse longer than 30 minutes.

§ Calibration or coding not required by the user.

|| Based on the sensor being replaced once every 14 days, and scanned at least once every 8 hours.

¶ The reader can capture data from the sensor when it is within 1 cm to 4 cm of the sensor.

# For Medicare patients, please visit <https://www.freestyle.abbott/us-en/medicare.html> for more information.

\*\*Participating pharmacies are subject to change without notice. Product availability may vary by retailer.

**References:** 1. FreeStyle Libre 14 day User's Manual. 2. Haak, Thomas, et al. "Flash Glucose-Sensing Technology as a Replacement for Blood Glucose Monitoring for the Management of Insulin-treated Type 2 Diabetes: a Multicentre, Open-label Randomised Controlled Trial." *Diabetes Therapy* 8, no. 1 (February 2017): 55-73. <https://doi.org/10.1007/s13300-016-0223-6>.

#### Important Safety Information

Failure to use FreeStyle Libre 14 day system as instructed in labeling may result in missing a severe low or high glucose event and/or making a treatment decision, resulting in injury. If readings do not match symptoms or expectations, use a fingerstick value from a blood glucose meter for treatment decisions. Seek medical attention when appropriate or contact Abbott at 855-632-8658 or [FreeStyleLibre.us](http://FreeStyleLibre.us) for safety info.

The product images are for illustrative purposes only.

The circular shape of the sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott.

© 2021 Abbott. ADC-45793 v1.0 10/21